

CELLiST™ BASAL Media (Model No. BASAL3)

Properties

- Chemically-defined, protein-free medium using no animal components.
- Contains no hydrolysates, extracts or other undefined components.
- Does NOT contain thymidine and hypoxanthine.
- Does NOT contain L-glutamine sources.
- Does NOT contain sodium bicarbonate and poloxamer.
- Contains 10.0 g/L D-glucose.

Instructions for use (preparation)

(1,000 mL of medium preparation)

1. Add all amount of BASAL3 powder (27.0g) to 900 mL of cell culture grade water (room temperature).
2. Put small amount of cell culture grade water in the pouch to rinse, add the water to 1 to ensure all the BASAL3 powder is used.
3. Quickly diffuse the powder into the water and mix until dissolved (approximately 30 minutes).
4. Add 1,000 mg of poloxamer P188 and 1,800 mg of sodium bicarbonate.
5. If necessary, add L-glutamine source to the solution at a proper concentration.
6. Mix until dissolved completely (approximately 30 minutes).
7. Adjust medium volume to 1,000 mL with cell culture grade water and mix to uniform consistency.
8. Sterilize by filtering through a membrane filter with pore size of 0.2 to 0.22 µm in diameter.
9. Keep prepared medium refrigerated (2 to 8°C) until use.
10. If necessary, add appropriate sterile growth factors (e.g. insulin or IGF-1) before use.

Storage

- Store under a cool (2 to 8°C), dark and dry condition.

FOR LABORATORY USE ONLY
NOT INTENDED FOR HUMAN OR THERAPEUTIC USE.

AJINOMOTO Co., INC.

Amino Acid Department, AminoScience Division
15-1, Kyobashi, 1-chome, Chuo-ku, Tokyo 104-0031, JAPAN
TEL: +81-3-5250-8111